

NEW PRODUCT RELEASE !



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SHOWER CHAIR

RG5301

FEATURES

FLAT PACK - Excellent for transport + storage
Easy to assemble - NO tools required

Contoured Plastic seat with draining holes
Height adjustable to suit every user
Rubber Leg tips for excellent stability in wet areas

Strong lightweight construction - Aluminium frame
Epoxy Powder coated framework
Easy clean design

Lifetime warranty against corrosion

SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions (at lowest setting):
565mm (W) x 460mm (D) x 770mm (H)

Seat Width:
480mm (between armrests)

Seat Depth:
420mm

Seat Height:
430 - 555mm (from floor to top of seat)
25mm adjustment increments

Armrest Height:
180mm (from top of seat to top of armrest)

Product Weight:
3kg

Maximum User Weight:
136kg

Product Materials:
Seat: Plastic (Polypropylene)
Frame / Legs: Powder coated Aluminium
Leg Tips: Flexible PVC

Minimum Order Quantity (MOQ):
4 units per box

Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

F: 08 9248 9458

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au