

# Dual Layer Ring Cushion



Unique central cavity for pressure free sitting

Designed to relieve pain in the tailbone and coccyx area

Standard cover

Waterproof cover

Dual layer for added comfort



A deluxe ring cushion with a dual layer construction for extra comfort.

Designed to relieve pain in the tailbone and coccyx area, this ring cushion comfortably conforms to body contours to reduce pressure point discomfort. Made of double layer, puncture-proof latex foam rubber, the affected areas are suspended by the unique central cavity for pressure free comfort when sitting.

Ideal for post natal or post operative discomfort, haemorrhoids, hip bursitis, coccyx or prostrate inflammation; or other painful conditions.

Available with either a standard stretch cover or waterproof cover. Machine washable at 30°C, the lightweight stretch cover offers added comfort, hygiene and ease of cleaning.

Suited to the incontinent user, the waterproof Polyurethane (PU) fabric cover also helps reduce sliding, friction and shear.



COVER	Standard	Waterproof
CODE	BL9566	BL9559
HEIGHT	80 mm	80 mm
WIDTH	400 mm	400 mm
DEPTH	340 mm	340 mm
CAPACITY	100 kg	100 kg
WEIGHT	0.5 kg	0.56 kg
WARRANTY	1 year	