

## Silicone Fibre Overlay



Cotton/polyester breathable cover provides maximum cool comfort



Silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention

This unique overlay helps to protect all pressure-prone areas of the body by reducing pressure, friction and shear, while helping to make any bed more comfortable and promoting good sleep.

The overlay is constructed from silicone coated cylindrical fibres that retain their loft and promote airflow to the skin's surface for superior patient support and comfort. The silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention, making it an ideal choice for geriatric or palliative care.

Covered with a breathable cotton/polyester fabric, the overlay is fully machine washable with standard detergent.

- fluid resistant silicone fibres provides extra support and comfort
- ideal for geriatric and palliative care
- protects pressure prone areas of the body
- flame retardant 16 CFR Part 1632
- fully machine washable up to 70°

<b>CODE</b>	BL9962
<b>HEIGHT</b>	100 mm
<b>WIDTH</b>	910 mm
<b>DEPTH</b>	1900 mm
<b>CAPACITY</b>	100 kg
<b>CONFORMS TO</b>	ARTG176835
<b>WARRANTY</b>	1 year

