

# MAVERICK 10




**ScooterWorld**  
AUSTRALIA



**merits**

*with you for life*



[www.scooterworldaustralia.com.au](http://www.scooterworldaustralia.com.au)

# MAVERICK 10



## FEATURES & BENEFITS

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Low centre of gravity design delivers a smooth and stable ride
- Swing away joy stick
- Height and width adjustable armrests
- Height adjustable seat with reclining back rest and swivel

Speed	8 kph (programmable)
Range	20 - 30km
Motor	DC 24 Volt / 520W
Controller	Dynamic Shark 40A or PG-VR2
Battery	12V/30-40 a/h x2 pieces
Charger	5 amp off-board
Gradient	6°
Front/Rear Wheel	6" x 2" foam filled
Drive Wheel	10" foam filled
Weight capacity	136kg
Brake	Intelligent, regenerative, electromagnetic brake
Length (Footplate down)	98 cm
Width	61 cm
Seat Width	46 cm
Seat Height (from deck)	43-53 cm
Seat Height (from ground)	51-56 cm
Turning Radius	50cm