

MAVERICK 14



ScooterWorld
AUSTRALIA

merits

with you for life

www.scooterworldaustralia.com.au



MAVERICK 14



FEATURES & BENEFITS

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension that works
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Low centre of gravity design delivers a smooth and stable ride
- Swing away joy stick
- Height and width adjustable armrests
- Height adjustable seat with reclining back rest

Range	25-35km
Motor	DC 24 Volt / 575W
Controller	Dynamic Shark 50A
Battery	12V/40-55 a/h x2 pieces
Charger	24V/5 amp off-board
Gradient	8°
Fron/rear wheel	8" x 2" foam filled
Drive Wheel	14" foam filled
Weight capacity	205kg
Brake	Intelligent, regenerative, electromagnetic brake
Length (Footplate down)	101 cm
Width	66 cm
Seat Width	51 cm
Seat Height (from deck)	46-51 cm
Seat Height (from ground)	56-58 cm
Turning Radius	53cm