

Introducing CairMax by Novis
– a non-powered hybrid foam air
combination pressure redistribution surface. It
is suitable for the prevention and treatment of skin
breakdown and pressure injuries in patients at risk.

			Standard	King Single
	CODE		FAMCM-R02	FAMCM-R02K
General	CAPACITY (SAFE WORKING LOAD)		230 kg	
	WARRANTY		2 year	
	PRESSURE INJURY RISK		At risk to high risk	
	ARTG		292670	
	CARTON DIMENSIONS (L x W x H) in mm		2270 x 900 x 180 mm	2270 x 1070 x 180 mm
Support Surface	DIMENSIONS (L x W x H) in mm		2000 x 850 x 200 mm	2000 x 1050 x 200 mm
	MATTRESS WEIGHT		15 kg	20.4 kg
		NUMBER	10	
		FOOT ZONE	Tapered Foam	
	CELLS	CELL MATERIAL	TPU laminated nylon	
		CELL FOAM INSERT MATERIAL	High resilience PU	
		MATTRESS AIR SUPPLY LINE	Heavy duty PVC piping	
	TOP COVER	MATERIAL / TYPE	PU-laminated nylon	
		FIRE RETARDANCY	BS7175 Crib 5	
		MAXIMUM WASH TEMPERATURE	95° C	
		ATTACHMENT	Dual zippers	
		SEAM CONSTRUCTION	External seams RF welded	
	BASE			

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.

