

QUEENSBURY



THEOREM CONCEPTS

RECLINER AND LIFT CHAIR SPECIALISTS

Shown in Santa Barbara
Latte Fabric

WEIGHT CAPACITY	158 kgs (350 lbs.)
SEAT-TO-FLOOR	46.9cm (18.5")
SEAT DEPTH	53.3 cm (21")
SEAT WIDTH	58.4 cm (23")
CHAIR DEPTH	105.5 cm (41.5")
CHAIR WIDTH	80 cm (31.5")
CHAIR HEIGHT	104 cm (41")
OPERATION	USB Hand Wand

The Queensbury delivers comfort and relaxation that is fully customisable to meet your needs. This contemporary four-motor power lift recliner offers an exceptional number of power-adjustable positions, including headrest, lumbar, backrest, extended footrest, and lift. A convenient remote with soft-touch buttons regulates each function independently. Available in two elegant fabric options, the Queensbury complements any living space.

FABRIC OPTIONS



SANTA
BARBARA
LATTE



PICASSO
WALNUT

FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  FOOTREST EXTENSION
-  POWER BACKREST
-  USB REMOTE

Viva LIFT!
POWER RECLINERS
By PHS® Mobility Products Corp.

Keep your Independence with...
SCOOTERS
& Mobility