

premium

CONCAVE POSITIONING COVER

- » a valuable accessory for the Premium 8, Premium 5 or Premium Slimline 3
- » helps prevent patients falling from the bed or injuring themselves from rolling against side rails
- » modular foam sections easy to remove for cleaning and patient transfer, while the hinged sections conform to profiling beds with head raises and knee breaks
- » components are heat sealed to protect against moisture ingress for maximum infection control.



| | | |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------|
| Code | TC-F358-PA | Concave Positioning Cover |
| | PTAM0051 | Concave Positioning Cover with Premium 8 |
| | PTAM0031 | Concave Positioning Cover with Premium 5 |
| | PTAM0011 | Concave Positioning Cover with Premium Slimline 3 |
| Warranty | 2 year | |
| Dimensions | L 2000 W 800 D 150 mm | |
| Approvals |     | |

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.

INFORMATION CORRECT AT TIME OF PUBLICATION - FEBRUARY 2017 © NOVIS HEALTHCARE, ALL RIGHTS RESERVED