

NEW PRODUCT RELEASE!











IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SOFT BRAKE **SEAT WALKER**

RG4204 - SILVER RG4204 - BLUE

FEATURES

Strong lightweight construction
Aluminium frame
User friendly Loop Brakes with soft brake lock
Simple brake cable adjustment
Height adjustable Handles
Ergonomic Handgrips to reduce pressure
Comfortable Padded Seat
Padded Backrest
Storage Bag
Large diameter wheels for excellent mobility
Flat free rubber tyres
Swivelling front wheels
Folds up for transit or storage

SPECIFICATIONS

Overall Dimensions (At Lowest Setting): 580mm (W) x 820mm (H) 635mm (Depth - Unfolded), 220mm (Folded)

Handle Height (From Floor to top of Handle): 820 - 920mm (5 settings)

Seat:

440mm (W - Between Handles) x 320mm (D) 540mm (From floor to Top of Seat)

Wheels

8" / Ø200mm (Front - Swivel / Rear - Fixed)

Product Weight:

7.5kg

Maximum User Weight:

136kg

Product Materials:

Contact Redgum brand for details

Please contact **REDGUM** if you require any more information on this product.

P: **08 9248 4180** F: **08 9248 9458**

E: orders@shoprider.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090 Part of the **FOR-DE Group** Pty Ltd.

redgumbrand.com.au