



NOW WITH INDICATORS

ADVENTURER

889SM EXCLUSIVE MODEL

NEW 2025



FEATURES

- Full adjustable suspension
- Comfortable high back contoured seat with swivel function
- Height adjust seat with position adjustment away from tiller
- Width and angle adjustable flip up padded armrests
- Height adjustable headrest
- Angle adjustable backrest
- Variable high/low speed settings
- Battery power indicator
- USB-A and USB-C Ports
- Gas assist adjustable tiller
- LED Headlight
- Rear lights
- Rear wheel drive
- Large pneumatic tyres
- Left and right indicators

SPECIFICATIONS

- Overall (lowest settings):
1300mm (L) x 640mm (W) x 1200mm (H)
- Seat:
520mm (W) x 520mm (D) x 620mm - 680mm (H)
(3 settings at 25mm increments)
- Backrest: **460mm (W) x 500mm (H) x 0 - 100°**
- Headrest:
255mm (W) x 180mm (H) x 630 - 750mm (H from seat)
(7 settings at 20mm increments)
- Armrests: **300mm (L) x 85mm (W) x 0 - 90° (Angle range)**
- Tiller: **460mm (W) x 900mm (H)**
x 30° forwards to 15° towards seat
- Turning Circle: **1500mm**
- Speed: **10km/h**
- Maximum range: **Up to 40km***
- Safe climb angle: **9°** / Maximum climb angle: **9°**
- Motor: **1.5Hp high torque motor (4 pole)**
- Batteries: **2 x 12V / 50 Amp (supplied)**
- Charger: **8 Amp (Supplied)**
- Tyres (Pneumatic): **Front / Rear: Ø330 x 100mm (4.00-5)**
- Product Weight: **99kg (with batteries)**
- Maximum User Weight: **159kg (on flat ground)**

* dependant on occupant and conditions.

RED
889SMR

SILVER
889SMS

ORANGE
889SMO

*Incline
weight rating*



0° - 159kg
6° - 113kg
9° - 91kg



ASK FOR THE BRAND

shoprider.com.au